



# ATOMIC SPECTROSCOPY

CEB 4032: ANALYTICAL CHEMISTRY

CFB3032: ANALYTICAL INSTRUMENTATION

TS. CHM. DR. MOHD DZUL HAKIM WIRZAL

Chemical  
Engineering

Inspiring Potential • Generating Futures

# Recap: Topics under Spectroscopy

## ***Atomic Spectroscopy***

- Atomic Absorption
- Atomic Emission
- Atomic Fluorescence
- Use of X-rays

## ***Molecular Spectroscopy***

- Fourier Transform Infra-red (FT-IR)
- Organic UV-Vis Absorption
- Effect of Solvent
- Inorganic UV-Vis Absorption
- Qualitative Analysis
- Quantitative Analysis
- Spectrophotometric Titrations
- Molecular Fluorescence
- Equipment for Molecular Spectroscopy

# At the end of the topic:

---



- 1) **Instrumentation** of atomic spectroscopy.
- 2) **Principle operation, differences, advantages and disadvantages** between AAS, AES and AFS.
- 3) **Quantitative analysis** of an **element** in an unknown sample.







# Introduction

- In **atomic spectroscopy**, samples vaporized at **2000-6000 K** and the **atomic concentration** are determined by measuring **adsorption, emission or fluorescence** at characteristic wavelength.
- **High sensitivity** → ability to distinguish one element from another in a complex sample.
- Ability to perform simultaneous **multi-element analyses**.
- Measuring **analyte concentration** at the parts-per-million level, parts per trillion.
- Precision → 1-2%, expensive, widely available.

# Overview: Atomic Emission, Atomic Absorption and Atomic Fluorescence

## Types of atomic spectroscopy:





(a), (b) and (c) represent the simplest possible atomic transitions that are used for each technique.

(d), (e), (f) represent the basic experiment arrangement that allows the transition to be observed

Atom cell: device that provides a **sufficiently high temperature** to break down the sample into its constituent atom → flame

# Outline

- Atomic Absorption
- Atomic Emission
- Atomic Fluorescence
- Use of X-rays



# Atomic Absorption Spectroscopy (AAS)


- Quantitative determination of chemical elements employing the **absorption of radiation (light)** by free atoms in the gaseous state.
- Used for determining the **concentration of a particular element (the analyte)** in a sample.



# AAS - Principle Operation

- Solution is aspirated into a flame and the sample element is converted to **atomic vapor**.
- Flame, contains atoms of that elements, some are thermally **excited by the flame**, but most remain in the **ground state**.
- Ground-state atom → **absorb radiation** of a particular wavelength .


# Principle Operation of AAS



- The absorption follows Beer's law:


$$A = k \cdot C$$

- A is absorbance, C is concentration of atomic vapor in the flame and **k is a specific coefficient** to each element at the given wavelength. The molar coefficient ( $\epsilon$ ) of absorption is not calculated in this case.
- Concentration of atomic vapor is directly proportional to the concentration of the analyte in the solution.
- Calibration curve of **absorbance** versus **concentration** in the standard solution → determination of concentration of the unknown solution.




Example of calibration graph in AAS


# AAS-Instrumentation



- It contains **five principal components**:
  - (1) light source
  - (2) burner/thermal device (element is brought to its atomic state)
  - (3) monochromator
  - (4) detector
  - (5) Readout device

# (1) Light source: Hollow-cathode Lamp (HCL)





Hollow-cathode lamps from an  
atomic absorption spectrometer






# (1) Light source: Hollow-cathode Lamp (HCL)

- Filler gas is ionized → positive ions are accelerated toward the cathode → **sputter metal atoms** from the cathode into gas phase.
- The free atom are **excited by collisions** with high-energy electrons and then **emit photons** to return to the ground state.
- The emitted energy can passed through the flame and can become **absorbed by the test element** because they posses the right energy.
- Hollow cylinder cathode coated with **element same as analyte**.
- Drawback: different lamp is required for each element.


# Principle of Ionization in HCL



## (2) Atomization: Burner / Thermal Device

### (i) Burner(Flame)

(a)




- Use a **premix burner** → fuel, oxidant and sample are mixed before introduce into flame.
- Solutions are aspirated with the gas flow into a nebulising/mixing chamber to form **small droplets** prior to entering the flame.
- Degree of atomization → temperature dependent


Fuel	Oxidant	Temperature (K)
Acetylene	Air	2400-2700
Acetylene	Nitrous oxide	2900-3100
Acetylene	Oxygen	3300-3400
Hydrogen	Air	2300-2400
Hydrogen	Oxygen	2800-3000
Cyanogen	Oxygen	4800

## (ii) Furnace

- Graphite furnace comprising **a tube of graphite** with a small cavity that can hold a precise quantity of sample  
→ requires less sample.
- Can attain **3000 K**.
- Graphite tube → surrounded by **double sleeve**
  - 1) inert gas (protect elements from oxidation)
  - 2) coolant (water).
- **Flameless** atomization → a very **high atom density** and a longer confinement period → **greater sensitivity of 1000**.



Solid sample weighed onto graphite platform




Direct solid sampling  
— end view of furnace



## (3) Monochromator

- Enables the **elimination of a large part of the stray light** due to the hollow cathode lamp, and the **selection of the most intense spectral line** in order to obtain a **better sensitivity**.






# Interference in Atomic Spectroscopy

- Any effect that **changes the signal** while **analyte concentration remain unchanged**.
- Type of interferences:
  - (1) **Spectral interference**
    - unwanted signals overlapping analyte solution.
  - (2) **Chemical interference**
    - chemical reactions decreasing the concentration of analyte atoms.
  - (3) **Ionization interference**
    - ionization of analyte atoms decreases the concentration of neutral atoms

# Background Correction

- To distinguish analyte signal from absorption, emission, and optical scattering of the sample matrix, the flame, plasma or furnace.
- To avoid significant error.
- Instrument with burner → low background noise.
- Most critical for graphite furnace (incomplete atomization).




Example of graphite  
furnace absorption  
spectrum


# Correction Methods

## (1) Beam Chopping (correct beam for flame emission)


- HCL (pulsing it on and off) → distinguish the signal of the flame from the **desired atomic line** at the same wavelength.
- the lamp being **periodically blocked** by the rotating chopper.
- while the **beam is blocked**, signal → from the **flame emission**.
- when the **beam is not blocked**, signal → from the **lamp and the flame**.
- **difference** between these signals → **desired signal**.



(a) Rotating chopper



(b)



(c)

**Figure 22-19** Operation of a beam chopper for subtracting the signal due to flame background emission. (a) Lamp and flame emission reach detector. (b) Only flame emission reaches detector. (c) Resulting square wave signal.


## (2) Deuterium Lamp Background Correction

- Additional  $D_2$  lamp passed through the flame in alternation with HCL.
- Light from HCL is absorbed **analyte and absorbed scattered** by background.
- Light from  $D_2$  lamp is absorbed and scattered **only by background**.
- **Difference** between absorbance measured  $\rightarrow$  **absorbance due to analyte**.


Disadvantage: uses more complex optical assembly including a second source



- 1 - Hollow cathode lamp
- 2 - Deuterium lamp
- 3 - Mirror
- 4 - Beam chopper (rotating mirror assembly)
- 5 - Grating
- 6 - Beam splitter
- 7 - Quartz window




With HCL



absorption:  
element and backg.

With D<sub>2</sub> lamp




absorption:  
backg. only

AAS with deuterium lamp background correction (double-beam AAS)

### (3) Smith-Hieftje (pulsed HCL) Correction

- When HCL run at high current  $\rightarrow$  output is broadened.
- The lamp is first run at low current, and pulsed at high current to measure absorbance.
- **Difference** between absorbances measured before and after  $\rightarrow$  **absorbance due to analyte.**




Disadvantage: high-current pulses reduce lamp lifetime




## (4) Zeeman Effect

- Shifting of energy levels of atoms and molecules in a **magnetic field**.
- **Sample and background** are observed when magnetic field is off.
- **Background alone** is observed when **magnetic field** is on.
- The difference → signal

Disadvantage: expensive



# Modern AAS with background correction



# Outline

- Atomic Absorption
- **Atomic Emission**
- Atomic Fluorescence
- Use of X-rays



# Atomic Emission Spectroscopy (AES)

- method of **chemical analysis** that uses the intensity of **light emitted** from **plasma, arc, or spark** at a particular wavelength to determine the **quantity of an element** in a sample.

## Advantages:

- Lower susceptibility** to chemical interferences.
- Good emission spectra → **multi-element analysis** for small sample.
- Energetic plasma source → determination of **low concentration of element** and **nonmetal** such as Cl, Br, I and S.


# Recall:





# AES-Principal

## Atomic Emission Spectroscopy





## Disadvantages:

- Emission spectral → made up hundreds or even thousands of lines → increase the probability of **spectral interferences** in quantitative analysis.
- Require **higher resolution and more expensive** optical equipment that needed for AAS.

## AAS

- simpler and less expensive
- lower operating cost
- greater precision
- require less operator

# AES- based on Arc and Spark


- Based on **excitation of emission spectra** of elements with **electric arcs or high voltage sparks** (Temp = 3000 - 6000 K).

## General operation:


- An electric arc or spark is passed through the sample → heating it to a high temperature **to excited the atoms** within it.
- The excited analyte atoms **emit light** at characteristic wavelengths that can be dispersed with a monochromator and detected.

# AES-based on Plasma Source

- Most common – Inductively Coupled Plasma (ICP-AES)
- Plasma- is a gas that is hot enough to contain **ions and free electrons**.
- Temperature up to **8000 K** → high enough to excite most elements



**A typical ICP source**



**Temperature profile of ICP**



**Photo of argon plasma in operation**







# Modern ICP-AES



# Outline

- Atomic Absorption
- Atomic Emission
- **Atomic Fluorescence**
- Use of X-rays

# Atomic Fluorescence Spectroscopy (AFS)



# Atomic Fluorescence Spectroscopy (AFS)

- can be considered a combination of AAS and AES in that the excitation from the ground state to the upper state is carried out radiationally (as in AAS), but the measurement is made by detection of radiation emitted by the atom as it relaxes back from the excited state to the ground state (as in AES).
- In AES, the atoms are excited only by the collisional energy inherent in the high temperature of the atom cell.
- In AAS and AFS, the atoms are excited primarily by radiation from a light source.

## Advantages:

- **Lower detection limit** → small quantity of fluorescence light can be more easily measured.
- **Greater sensitivity** achievable because the fluorescence signal has a very low background.

## Disadvantages:

- **High cost and operational complexity** in AFS; and **great success and widespread application** in AAS → discouraged further development in AFS.



## Modern AFS





# Detection Limit for Atomic Spectroscopy

Element	AAS Flame	AAS Electro- thermal	AES Flame	AES ICP	AFS Flame
Al	30	0.1	5	0.2	5
As	200	0.5	—	2	15
Ca	1	0.25	0.1	0.0001	0.4
Cd	1	0.01	2000	0.07	0.1
Cr	4	0.03	5	0.08	0.6
Cu	2	0.05	10	0.04	0.2
Fe	6	0.25	50	0.09	0.3
Hg	500	5	—	—	5
Mg	0.2	0.002	5	0.003	0.3
Mn	2	0.01	—	0.01	1
Mo	5	0.5	100	0.2	8
Na	0.2	0.02	0.1	0.1	0.3
Ni	3	0.5	600	0.2	0.4
Pb	8	0.1	200	1	5
Sn	15	5	300	—	200
V	25	1	200	0.06	25
Zn	1	0.005	50000	0.1	0.1



# Application of Atomic Spectrometry

- Widely used in many laboratories → when trace **element analyses** are required.
- Environmental samples for **heavy –metal** contamination.
- Pharmaceutical samples → analyzed for **metal impurities**.
- Steel industry → determine **minor and major components**.

Choice of atomic spectrometry method → Depend on the sensitivity required, number of samples, **single-element or multi-element measurements**.

# Example 1

Five standard solutions were prepared for measuring the lead concentration in two solutions, A and B. The two solutions **A and B contain the same concentration of magnesium** used as an **internal standard**. The following data were obtained.


Conc. (mg/L)	Signal of Pb	Signal of Mg
0.10	13.86	11.88
0.20	23.49	11.76
0.30	33.81	12.24
0.40	44.50	12.00
0.50	53.63	12.12
Sol A	15.50	11.80
Sol B	42.60	12.40

From the data, calculate the **lead concentration (mg/L)** in the two sample solutions, A and B.

# Solution

Construct the **calibration graph** of the **ratio of the signal emission Pb/Mg** against **concentration**

Conc. (mg/L)	Signal of Pb	Signal of Mg	Ratio Pb/Mg
0.10	13.86	11.88	1.17
0.20	23.49	11.76	2.00
0.30	33.81	12.24	2.76
0.40	44.50	12.00	3.71
0.50	53.63	12.12	4.42
Sol A	15.50	11.80	1.31
Sol B	42.60	12.40	3.44



Ratio Pb/Mg,  $y = 8.2274 x + 0.3437$ ,  
 whereby  $x$  = Pb concentration of solution in mg/L

For Solution A, Ratio **Pb/Mg = 1.31**

$$1.31 = 8.2274 x + 0.3437$$

$$x = 0.12 \text{ mg/L}$$

Therefore, concentration of **Pb in Solution A = 0.12 mg/L**

For Solution B, Ratio **Pb/Mg = 3.44**

$$3.44 = 8.2274 x + 0.3437$$

$$x = 0.38 \text{ mg/L}$$

Concentration of **Pb in Solution B = 0.38 mg/L**

## Example 2

A stock solution of calcium ions was prepared by dissolving **1.834 g of  $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$**  in **1.00 L** of distilled water and then further diluting by a **factor of ten**.

From this new solution, three standard solutions were prepared by further dilutions of **five, ten and twenty times**, respectively. The **unknown sample** is itself diluted **twenty-five times**.

Sufficient strontium chloride was then introduced to eliminate any interference due to phosphate ions. An analytical blank containing the same concentration of strontium was the first solution to be examined by the air/acetylene flame. The results were as follows:



Standard	Signal
Blank	1.50
1:20	10.6
1:10	20.1
1:5	38.5
Unknown	29.6

MW:

Ca = 40.1 g/mol

Cl = 35.5 g/mol


What is the **concentration (in ppm) of the calcium** in the unknown sample?

# Solution

- 1) The molar mass of  $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$   
 $= (40.1 + (35.5 \times 2) + (18 \times 2)) = 147.1 \text{ g}$
- 2) The concentration of the parent solution (in 1.00 L solution)  
 $= 1.834 / 147.1 = 0.01247 \text{ mol/L}$
- 3) Solution diluted 10 times  $= 0.01247 / 10 = 0.001247 \text{ mol/L}$
- 4) Tabulated the data

Standard	Signal	Concentration (mol/L)
Blank	1.50	0
1:20	10.6	$0.6234 \times 10^{-4}$
1:10	20.1	$1.2468 \times 10^{-4}$
1:5	38.5	$2.4935 \times 10^{-4}$
Unknown	29.6	?

## 5) Plot signal versus concentration



## 6) Concentration of unknown sample (after dilution)

Given signal of the unknown sample  $y = 29.6$

$$x = (29.6 - 1.46)/148635 = 1.8932 \times 10^{-4} \text{ mol/L}$$

7) The parent solution has a concentration **twenty-five times greater**, therefore the concentration of unknown solution

$$= 1.8932 \times 10^{-4} \text{ mol/L} \times 25 = 4.7331 \times 10^{-4} \text{ mol/L}$$

8) Concentration of  $\text{Ca}^{2+}$

$$= 4.7331 \times 10^{-4} \text{ mol/L} \times 40.1 \text{ g/mol} = 0.19 \text{ g/L}$$

$$\begin{aligned} \text{In ppm} &= 0.19 \text{ g/L} \times 1\text{L}/1000 \text{ mL} \times 10^6 \mu\text{g/g} \\ &= \underline{190 \mu\text{g/mL} = 190 \text{ ppm} = 1.90 \times 10^2 \text{ ppm}} \end{aligned}$$

# Example 3

An AAS method is employed for the determination of lead (Pb) in a sample of adulterated paprika by the introduction of lead oxide (of same colour). An electrothermal atomic absorption instrument that provides a background correction based upon the Zeeman effect is used.

**0.0100 g** of the paprika powder is placed in the tube of the graphite furnace. The determination of the area peak absorbance was made at  $\lambda = 283.3$  nm first in the absence and later in the presence of a magnetic field. The value of the peak absorption following background correction was **1220** (arbitrary units).

Under the same conditions, **0.0100 mL of a solution of 10.0 g/L Pb** led to a value of **1002** in the same units.

Calculate the **% mass of lead** in the sample of paprika under study.

# Solution

The signal of absorbance of the sample extracted from paprika =  $A_x$  and the sample solution =  $A_s$ ,

$m_x$  = mass of lead in the paprika

$m_s$  = mass of lead in the sample solution

We have:  $A_x/A_s = m_x/m_s$ , therefore  $m_x = m_s A_x/A_s$

Given:  $A_x = 1220$ ,  $A_s = 1002$ ,

$$m_s = 10.0 \text{ g/L} \times 0.0100 \text{ mL} \times (1\text{L}/1000 \text{ mL}) = 0.0001 \text{ g}$$

$$m_x = (0.0001 \text{ g} \times 1220/1002) = 0.000122 \text{ g}$$

The % mass in the paprika sample is therefore  
 $= (0.000122/0.01) \times 100\% = \underline{1.22\%}$



# Exercise 1

Magnesium is determined in a river water sample using electrothermal AAS, using multiple standard additions. Four 25.0 mL aliquots are pipetted into 50.0 mL volumetric flasks, and aliquots of a 2.50 ppm standard are added in the amounts of 0.00, 1.00, 2.00, 3.00 mL, respectively. Fixed volumes of each solution are measured, giving net absorbance signals of 0.10, 0.18, 0.24, 0.31, respectively. Calculate the concentration of magnesium in the sample.


# Outline

- Atomic Absorption
- Atomic Emission
- Atomic Fluorescence
- Uses of X-rays

# Atomic X-rays Spectroscopy

- X-rays are short-wavelength of electromagnetic radiation produced by the deceleration of high energy electrons or by electronic transitions in the inner orbitals of atoms.
- Based on measurement of emission, absorption, fluorescence and diffraction of electromagnetic radiation.
- The wavelength range:  $10^{-5}$  Å to 100 Å. Conventional X-ray spectroscopy: 0.1 Å to 25 Å
- X-ray methods are widely used for the qualitative and quantitative determination of elements in the periodic table.

# Recall



# X-ray

- X-rays are generated:

- (1) By **bombardment** of a metal target with a beam of high-energy electron
- (2) By exposure of a substance to a primary beam of X-rays to generate a **secondary beam of X-ray** fluorescence.
- (3) By use of a **radioactive source** whose decay process results in X-ray emission.

# X-rays Instrumentation

- Generally consists of **five main components**:
  - (1) **source** → tubes, radioisotopes and secondary fluorescence.
  - (2) **monochromator (collimators)** → restricting the wavelength range of incident radiation.
  - (3) **sample holder**
  - (4) **radiation detector/transducer** → convert radiant energy to electrical signal.
  - (5) **signal processor and readout**



# X-ray Fluorescence (XRF)

- The elements in the sample **excited by absorption of the primary beam** and **emit their own characteristic fluorescence X-rays**.

# Use of X-ray Fluorescence (XRF)


- Quality control in the **manufacture of metal and alloys**.
- **Correct the composition** of the alloy during its manufacture.
- Easily adapted to liquid sample → **direct quantitative determination** of lead and bromine in aviation fuel sample.
- Determination of Ca, Ba and Zn in the **lubricating oils**.
- Determination of **pigment paint**.
- Analysis of **pollutant**.
- Determination of elements heavier than Na in rocks and soil.

# Use of X-ray Absorption

- Widely-used technique for determining the **local geometric and/or electronic structure** of matter.
- Determination of **lead in gasoline** and **sulfur or halogen in hydrocarbons**.
- General useful only when a **single element** with high atomic number.
- Relatively **limited** when compared with XRF procedure.
- Cumbersome and time-consuming.


# Use of X-ray Diffraction (XRD)

- Important information to science and industry.
- Determine the **arrangement and the spacing of atoms** in crystalline materials.
- Understanding of **physical properties** of metals, polymeric materials, and other solids.
- Determine the **structures of such complex natural products** such as steroids, vitamins and antibiotics.
- **Qualitative identification of crystalline compounds** → each crystalline substance has a unique X-ray diffraction.
- X-ray powder diffraction → providing **qualitative and quantitative information** about the compounds present in a solid sample



**X-ray Diffractometer**

# X-ray Diffraction Spectrum





X-ray source: XR50

X-ray monochromator





X-ray Fluorescence Spectrometer



X-ray Absorption Spectrometer



X-ray Emission Spectrometer



## X-ray Diffraction Spectrometer





# END OF ATOMIC SPECTROSCOPY